

Rest Recovery PEMF Therapy System — Specification Sheet

Rest Recovery Studios LLC

Info@RestRecoveryWellness.com

(All specifications subject to change. For installation or permitting use only.)

Product Overview

The **Rest Recovery PEMF Therapy System** is a high-intensity, dual-program electromagnetic therapy device designed for both **clinical** and **home applications**. With intelligent software and a large touchscreen interface, it delivers adjustable pulsed or continuous fields for optimal performance across various treatment needs, including musculoskeletal recovery, pain relief, and circulation enhancement.

General Specifications

Model	Rest Recovery PEMF Pro Max
Application	Human and animal therapy
Power Consumption	850 W
Magnetic Intensity	Up to 24,700 Gauss
Pulse Frequency	1–10 Hz
Oscillation	4.5 kHz
Voltage	110 V AC
Amperage	15 A
IP Rating	IP31
Weight	35 kg
Dimensions (W x D x H)	75 x 58 x 52 cm
Display	12.1" Intelligent Touchscreen
Program Modes	Pulsed / Continuous
Preset Protocols	Clinical and veterinary presets included

Accessories

Dual Paddle	Thigh, back, hip, abdomen — Up to 19,800 Gauss
Butterfly Loop	Head, shoulders, knees, ankles, joints — Up to 24,700 Gauss
Large Loop	Chest, back, stomach, thighs — Up to 8,900 Gauss
Mattress	Full-body treatment — Up to 8,300 Gauss
Cube	Pelvic floor, hip pain, urinary incontinence — Up to 9,000 Gauss
X-Wings (Equine)	Broad-area equine therapy — Up to 8,550 Gauss
Hoof Box (Equine)	Hoof recovery and inflammation — Up to 9,000 Gauss

Electrical Requirements

Voltage	110 V Single Phase
Amperage	15 A Dedicated Circuit
Max Power Draw	850 Watts
Plug Type	NEMA 5-15P Grounded Plug
Breaker	Dedicated 15 A GFCI Breaker Recommended

Features

- Dual treatment programs (Human & Veterinary modes)
- Pulsed or continuous frequency operation
- Intelligent preset therapy protocols
- 12.1" smart touchscreen interface
- Adjustable magnetic intensity and pulse rate
- High output up to 24,700 Gauss
- Multiple applicator attachments for targeted or full-body use

Certifications

CE / ETL / ROHS / ISO 9001 Certified Components

Maintenance

- Wipe down handles and loops daily
- Inspect power cord and applicator connections weekly
- Ensure adequate ventilation around control unit
- Avoid excessive moisture or humidity

Warranty

1-Year Limited Warranty on console and applicators

(Warranty void if used on ungrounded power or outside specified amperage.)